OSTEOPLUS®

Aditamentos, componentes, instrumental e implantes dentales OSTEOPLUS® con Mount

Implantes dentales OSTEOPLUS®: Iridium, Argentum y Platinum.

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INFORMACIÓN IMPORTANTE

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OSTEOPLUS®

OSTEOPLUS® abutments, components, instruments and dental implants with Mount.

OSTEOPLUS® dental implants: Iridium, Argentum and Platinum.

INSTRUCTIONS FOR USE

IMPORTANT INFORMATION

Read these instructions for use carefully before using the product.

This document contains basic information for the use of the prosthetic instruments and abutments or components of OSTEOPLUS®, original dental implant systems, from now on referred to as OSTEOPLUS®. This document has been written as a quick reference guide for the physician in charge of the treatment, henceforth, “User”. It does not provide sufficient indications and technical specifications for the use of OSTEOPLUS® products. It is neither an alternative nor a substitute for specialised training and professional clinical experience.

OSTEOPLUS® products must be used carrying out sufficient planning of the treatment and rigorously following the surgical and/or prosthetic protocols set by Henry Schein España SL Read carefully the specific surgical and prosthetic protocol of each product, the Instructions of Use and the Maintenance Instructions before using a OSTEOPLUS® product.

1. Description of the system

OSTEOPLUS® products are composed of dental implants, prosthetic components or components and surgical and prosthetic instruments.

Henry Schein España SL has designed a wide range of dental implant systems for the carrying out of dentistry treatments in accordance with the existing scientific bibliography and current clinical standards. Henry Schein España SL commercialises latest-generation implant products and a wide range of restorative solutions to respond to the particular needs of each patient.

2. Description of the products

Dental implants

OSTEOPLUS® dental implants are manufactured in biocompatible grade 4 titanium, according to ASTM F67 and ISO 5832-3 standards, with Double Acid Engraving surface treatment. OSTEOPLUS® dental implants are manufactured in different diameters, lengths, shapes and connections depending on their intended purpose.

Iridium, Argentum and Platinum dental implants come with a metal mount, which makes it easier to transfer them to the patient’s mouth and start placement at the implant site.

Instruments

OSTEOPLUS® instruments are manufactured in stainless steel, grade 5 titanium and/or plastic. Some instruments are laser-engraved to facilitate identification and use.

OSTEOPLUS® instruments

Insertion keys- Screwdrivers- Screwdriver tips.

Abutments

OSTEOPLUS® prosthetic components or components are manufactured in grade 5 titanium, Cr-Co alloy and biocompatible plastic materials.

OSTEOPLUS® prosthodontic abutments or components

Castable straight abutments- Castable angled abutments- Transepithelial abutments - Abutments for overdentures - Interface abutments - Temporary abutments - Clinical screws- Healing caps- Healing abutments- Abutments with machined base - Impression abutments - Abutment adapters - Aesthetic abutments - Overdenture retention systems-- Laboratory screws- Analogues- Scanbodies- UCLA.

3. Indications for use

It is essential to previously evaluate the patient by means of a clinical and imaging diagnosis and study of models, as well as to determine the possible contraindications for the treatment, prior to carrying out a treatment with OSTEOPLUS® products. A correct treatment with dental implants must be based on a diagnosis, a treatment plan and a design of the prosthetic rehabilitation which meets the needs and expectations of the patient.

Implants

OSTEOPLUS® dental implants are indicated for both mandibular and maxillary rehabilitation and are intended to restore chewing function, as well as providing improved occlusion and distribution of occlusal loads, phonation, aesthetics, prosthetic stability and maintenance of the alveolar process, in addition to a positive psychological effect. They are supplied with different designs and in different lengths and diameters so that they can be adapted to different clinical cases.

OSTEOPLUS® dental implants are single-use products.

Important

Carefully consult the clinical indications for each implant model in the appropriate surgical catalogues and protocols before planning surgery. Do not use OSTEOPLUS® dental implants in dental positions or situations that are not described as indications in the surgical catalogues or protocols.

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Instruments

OSTEOPLUS® surgical and prosthetic instruments are designed as one of the fundamental parts of the implant treatment. With the instruments, you will be able to carry out all the necessary processes for the preparation of the bone tissue, the insertion of the OSTEOPLUS® dental implants and the manipulation of the different OSTEOPLUS® prosthetic components or components.

The surgical drills and trephine drills are designed to withstand a maximum of 45 uses under optimum conditions, whilst the circular scalpels are designed to withstand a maximum of 25. Improper use, lack of cleaning and maintenance or inadequate disinfection can reduce its number of uses.

OSTEOPLUS® instruments are reusable.

Abutments

OSTEOPLUS® prosthetic components or components are fixed to the OSTEOPLUS® dental implants and serve as a support for the prosthetic rehabilitation. The prosthetic rehabilitation may be single or multiple, partial or total, and may be screwed or cemented to the abutments, or may be removable by means of implant anchorages or retainers.

The attachments or abutments in biocompatible plastic material (PEEK) are for provisional use and are designed to be used in the mouth for a period of time no longer than 180 days. The removable prostheses with anchorages or retainers must be monitored at least once per year by the clinical professional. The correct intraoral adjustment of the prosthesis and the retention between the anchorages or retainers and the dental prosthesis must be verified. A loss of adjustment and/or a poor retention of the prosthesis put the prosthetic rehabilitation and the dental implants in danger.

OSTEOPLUS® abutments are single-use products.

4. Contraindications

It is necessary to carry out a preoperative medical examination of the patient and to determine if absolute and/or relative contraindications or risk factors for the treatment with dental implants exist. The user is responsible for evaluating the potential risks and benefits of the implant treatment and for taking the pertinent decisions depending on the case. OSTEOPLUS® products must not be used on patients who lack the minimum medical conditions for receiving a prosthetic treatment and rehabilitation with dental implants.

Relative contraindications

Age of the patient - Uncontrolled systemic disease - Poor oral hygiene - Active local infection- Pregnancy- Therapy with bisphosphonates- Sjogren syndrome- Uncontrolled endocrine disease - Radiation to head and neck - Allergy to stainless steel. - Stress - Psychological factors - Valve prostheses - Uncontrolled periodontal disease - Smoking - Alcoholism - Drug addiction - Addictions in general.

Absolute contraindications

The manufacturer does NOT recommend the implant treatment in the following cases: Cardiovascular diseases or disorders - Haemophilia - Cirrhosis of the liver - Terminal pathologies - Epilepsy - Neoplasms under care with chemotherapy- Blood dyscrasias - Cancer - Uncontrolled diabetes - Untreated periodontitis - Infectious diseases - Leukaemia - Alterations of the immune system- Treatment with immunosuppressives - Allergy to titanium - Active treatments with radiotherapy - Indefinite treatments with corticosteroids- Treatments with anti-coagulants- Involuntary muscle disorders - Serious liver disorders - Endocrine disorders - Psychiatric and/or psychological disorders.

The Manufacturer will update these contraindications, as necessary, based on information obtained from clinical studies and/or post-marketing information.

5. Storage

OSTEOPLUS® products must be stored in their original packaging, in a clean, dry place, protected from sunlight and adverse conditions. In the case of OSTEOPLUS® Implants they must also be preserved at a temperature between -10 and 55°C.

6. Packaging and supply state

Dental implants are supplied in a cardboard box with a product identification label, which acts as a tamper-evident seal. The carton has a shrink-wrapped plastic film that forms a safety seal to prevent handling.

This package contains:

- QR code to allow access to these instructions for use of the product.
- Heat-sealed plastic blister with label on the back, indicating the manufacturer, the complete reference of the product and the lot number. The implant mount (implant attached to metal mount) container is inside.
- Adhesive labels to allow implant traceability (these labels indicate the manufacturer, complete product reference number and lot number). One label should be adhered to the patient’s implant card while the others are for the clinic’s records.

OSTEOPLUS® instruments and prosthetic components or components are supplied packaged in a thermosealed plastic bag; the product inside is identified with a label stuck to the back.

Henry Schein España SL recommends that the traceability of the abutments used be maintained by recording the reference and the lot used in in the patient’s file.

Henry Schein España SL guarantees that all its products follow an exhaustive process of manufacture, control and cleaning before being packaged. The correct integrity of the packaging must be controlled before using the OSTEOPLUS® products, checking

that it is not damaged. If there is any defect or damage, it must not be used and must be notified to Henry Schein España SL.

7. Cleaning and sterilisation

OSTEOPLUS® dental implants are sterilised with radiation of Beta Rays at 25kGy. In the event that the sterile container is altered, handled, damaged or broken, sterility is not guaranteed and the product should not be used. Do not use the product after the expiry date indicated on the packaging.

The re-sterilisation of OSTEOPLUS® dental implants is NOT authorised.

The instruments and prosthetic attachments or components of OSTEOPLUS® dental implant systems are supplied NOT sterilised.

Instruments and prosthetic attachments or components that are to come into contact with the patient’s mouth must be cleaned, disinfected and sterilised prior to use in accordance with the manufacturer’s cleaning, disinfection and sterilisation instructions provided in the catalogues. Visit our website to view the latest version at www.henryschein.es

[For the US, the approved, validated sterilisation cycle uses a gravity-displacement steam autoclave at 132°C / 270°F for 15 minutes and a drying time of 15-30 minutes].

Important

Neither the inside nor the outside of the plastic blister pack containing the implant vial is sterile. Only the inside of the container is sterile; do not place the blister in the sterile surgical field. The sterile container should NOT be opened until just before implant placement.

8. Precautions

Surgery and Instruments

The surgical protocols describe in detail the precautions and important factors to take into account during the surgical process of fitting the dental implants. Prior to the intervention, it is essential to have consulted the surgical protocol of the product, as well as all the indications about its correct usage indicated in the catalogues of OSTEOPLUS®.

Particular attention must be paid and the necessary measures taken to prevent the instruments, of such small size, from being swallowed or breathed into an airway. All OSTEOPLUS® instruments for manual use incorporate retaining elements between each other or for use with dental floss, to avoid accidental ingestion.

It is important to reduce the time of surgery and the damage done to the area receiving the implants, paying special attention to bone trauma and the asepsis of the surgical field. The correct preparation of the surgical site is a fundamental factor in the success of the treatment. This process requires that the surgical protocol for drilling set by the manufacturer for each type of implant and the correct use of the surgical drills, both in terms of the number of uses and of the working speeds, be followed rigorously. The surgical drills and taps must be used with constant abundant irrigation to prevent overheating of the bone.

Take special care not to apply excessive forces or torques during the insertion of the implant. It is important to control the insertion torque of the implant, both in insertion with contra-angle, and with insertion with ratchet wrench.

Attention

When a surgical instrument with a limited number of uses is used, that number must not be passed, as this would put the process of osseointegration, and therefore the treatment, in grave danger.

The maximum torque for insertion of the dental implants is 50 Ncm. Exceeding the maximum insertion torque indicated for implants can cause serious damage to the dental implant in its connection, in the carrier and in the included clinical screw. Refer to the surgical protocol for specific carrier removal considerations based on implant connection type and bone type.

Prosthetic rehabilitation

The surgical protocols describe in detail the precautions and important factors to take into account during the surgical process of fitting the prosthetic rehabilitation. The design, the type of rehabilitation and the execution of the dental prosthesis must be studied prior to the insertion of the dental implants. Take special care not to apply excessive forces during the tests of the abutments or dental prostheses; this must always be done manually. Pay special attention to the definitive adjustment of the rehabilitation; the torque applied to the definitive clinical screws must be correctly controlled following the manufacturer’s indications, whether carried out mechanically or manually.

Table of OSTEOPLUS® abutment torques		
Element/Abutment	Instrument/Tool	Torque
Healing caps/healing abutments	Hex screwdriver 1.25 mm	Manual
Impression abutment screws	Hex screwdriver 1.25 mm	Manual
Laboratory screws	Hex screwdriver 1.25 mm	10 Ncm
Direct-to-implant clinical screws	Hex screwdriver 1.25 mm	30 Ncm
OsteoPlan® abutments	Drivers: xxxxxxx	30 Ncm
Clinical screws on OsteoPlan®	Hex screwdriver 1.25 mm	20 Ncm
OsteoCam® scanbody + screw	Hex screwdriver 1.25 mm	Manual
Osteolock® abutments	Drivers: 702993	30 Ncm

Attention

Applying insufficient torque can cause incidents of loosening of said screws. If this occurs, we recommend replacing them with new products, as retightening a screw increases the risk of fatigue and fracture of it.

Exceeding the maximum torque applied to the clinical screws indicated by the manufacturer can cause serious damage to the dental implant, to its connection, to the prosthetic abutment and to the clinical screw. Always check the passive adjustment of the implant-prosthesis unit. An incorrect passive adjustment puts the prosthetic rehabilitation in grave danger and can cause future damage to the dental implant, to its connection, to the prosthetic abutment, to the clinical screw and/or cause the loss or reabsorption of bone.

9. Warnings

The planning of the treatment and the fitting of OSTEOPLUS® dental implants requires specific dentistry training. Product specifications alone do not guarantee proper use of the product. The users are recommended to attend specialised theoretical/practical training courses to become familiar with the work techniques and protocols suitable for each product, including biomechanical, radiographical and prosthetic requirements associated with the implant treatment. It is the user's obligation to know and be informed of the current state of the technique for each product and its possible applications.

Before using OSTEOPLUS® products, it is necessary to know and be familiar with the corresponding surgical and prosthetic procedures. The patient must meet sufficient anatomical and psychological conditions to undergo the treatment with dental implants. The user must ensure the correct planning of the treatment, guaranteeing its safety with a minimum margin for error and respecting the vital oral structures and the general health of the patient.

The procedure for use of the surgical drills, taps and other surgical instruments needed for the fitting of the implants is specified in the surgical protocols corresponding to each type of implant. The fitting of the implant and the prosthetic planning must adapt to the patient's individual conditions, and in particular to a correct distribution of forces. The passive fit in the prosthetic rehabilitation structure and the adjustment of the occlusion must be achieved, and the appearance of excessive lateral forces avoided. A treatment with an insufficient number of implants, or an unsuitable choice of size and position to withstand and transmit the loads expected, can cause mechanical failure of the implant, of the abutment or prosthetic components or of the clinical screws due to overburdening or fatigue and substantial bone loss around the implant.

A qualitatively or quantitatively insufficient alveolar bone, the appearance of infections or diseases in general, and changes in the habits of the patient, are some potential causes of failure of the osseointegration and of the treatment. The lack of bone tissue or soft tissue can give rise to a deficient aesthetic result. An unsuitable prosthetic rehabilitation can cause the failure of the whole rehabilitation treatment.

The reuse of single-use products carries the danger of possible deterioration of their properties. The geometry of the product may have suffered damage during previous use or due to an unsuitable use. This implies the risk of surgical or prosthetic failure and/or harm to the patient's health.

Electrosurgery is not indicated due to the conductivity of the dental implants.

10. Adverse or side effects

Many adverse or undesired effects can present themselves in current implantology. These are documented in the specialised scientific bibliography and published in journals and books specialised in the dentistry sector.

However, the most relevant are:

Post-operative discomfort - Localised pain - Localised inflammation - Haematomas - Haemorrhages - Localised or systemic infections - Difficulty speaking - Bone loss and/or fracture - Loss of the implant - Damage to adjacent teeth - Fractures to the implants and/or prosthetic components - Fracture of instruments - Injury to adjacent anatomical structures (inferior alveolar nerve or maxillary sinus).

The Manufacturer will update these contraindications, as necessary, based on information obtained from clinical studies and/or post-marketing information.

11. Information about responsibility, safety and guarantee

The handling and use indication of OSTEOPLUS® products are based on published international literature, current clinical standards and our clinical experience with our products and should therefore be understood as indicative general information. The handling and use of OSTEOPLUS® products, as they are beyond the control of Henry Schein España SL, are the sole responsibility of the user. Henry Schein España SL, its subsidiaries and/or official authorised distributors assume no responsibility, whether express or implicit, full or partial, for the possible harm or damages which may result from the improper handling of the product or from any other event not specified in its protocols and manuals for the correct use of its products.

The User of the product must ensure that the OSTEOPLUS® product used, is suitable for the intended purpose and procedure. Neither these instructions for use nor the protocols of work or manipulation of the products exempt the user from the clinical obligation. The use, handling and clinical application of OSTEOPLUS® products must be carried out by qualified professional personnel with the necessary qualifications according to the current legislation of each country. The use, manipulation and/or application, in whole or in part, of OSTEOPLUS® products in any of their manufacturing phases by unqualified personnel or without the necessary qualifications, automatically voids any type of guarantee and may cause severe damage to the health of the patient.

OSTEOPLUS® products have their own system, with designed features and working protocols, including dental implants, prosthesis components and fittings, and surgical or prosthetic instruments. The use of OSTEOPLUS® products in combination with elements or components from other manufacturers can cause failure of the treatment, cause tissue damage, cause damage to the bone structures, inadequate aesthetic results and serious harm to the patient's health. For this reason, only original OSTEOPLUS® products should be used.

The clinician, who is responsible for the treatment, is solely responsible for ensuring that original OSTEOPLUS® products are used and that they are used in accordance with the corresponding instructions for use and catalogues throughout the entire implant treatment process. The use of non-original OSTEOPLUS® components, instruments or any other product used alone or in combination with any of the original OSTEOPLUS® products will automatically void any guarantee on the original OSTEOPLUS® products.

Warning

Non entire OSTEOPLUS® products are available in all countries. Please consult for availability.

The brand OSTEOPLUS® is a registered trademark, as some of the products mentioned or not in this document. Henry Schein España SL reserves the right to modify, change, eliminate or evolve any of the products, prices or technical specifications referenced in this document or any of its catalogues without prior notice. All rights reserved. The whole or partial reproduction or publication of this documentation, by any means and in any format, is forbidden without the corresponding authorisation in writing from Henry Schein España SL.

12. Magnetic resonance

The safety and compatibility of OSTEOPLUS® products in a magnetic resonance environment has not been evaluated, nor have they undergone testing to see if they heat up or migrate in said environment.

13. Removal of the material

The dental implants, abutments and kit, as well as each product and consumable used during surgery for placing dental implants and which are disposed of may pose a health hazard to people handling them. For this reason, you are advised to check and comply with current legislation on biohazardous waste.

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Symbols of labelling codes and used in catalogues, standard ISO 15223-1

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CE marking and notified body number

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A single sterile barrier system. Sterilised using irradiation

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A single sterile barrier system with a protective outer packaging. Sterilised using irradiation

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Caution, Consult accompanying documents

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OPENING OF THE PACKAGE AND HANDLING

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EN Open cardboard box by the flap and remove blister.

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EN Carefully remove the attached seal that seals the blister.

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
EN Hold the vial with one hand in a vertical position. Remove the cap by turning it vertically.

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EN Remember to remove the identification labels from the implant to stick them on the Implantation Document and on the patient history sheet in order to document the traceability of the reference and batch of the implant.

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